

MAP2K5 Mouse mAb

Catalog No: #63899

Package Size: #63899-1 50ul #63899-2 100ul


Description

Product Name	MAP2K5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MAP2K5 (AA: 63-180) expressed in E. Coli.
Target Name	MAP2K5
Other Names	MEK5; MAPKK5; PRKMK5; HsT17454
Accession No.	Q13163
Calculated MW	50kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using MAP2K5 mouse mAb against Jurkat (1), A431 (2), A549 (3) cell lysate.



Flow cytometric analysis of Jurkat cells using MAP2K5 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using MAP2K5 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.