

MAP2K6 Mouse mAb

Catalog No: #63901

Package Size: #63901-1 50ul #63901-2 100ul


Description

Product Name	MAP2K6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human MAP2K6 expressed in E. Coli.
Target Name	MAP2K6
Other Names	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; MAP2K6
Accession No.	P52564
Calculated MW	38kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using MAP2K6 mouse mAb against HepG2 (1), MCF-7 (2) and NIH/3T3 (3) cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer (left) and kidney cancer (right) using MAP2K6 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.