

MAP3K14 Mouse mAb

Catalog No: #63903

Package Size: #63903-1 50ul #63903-2 100ul


Description

Product Name	MAP3K14 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MAP3K14 (AA: 769-947) expressed in E. Coli.
Target Name	MAP3K14
Other Names	HS; NIK; HSNK; FTDCR1B
Accession No.	Q99558
Calculated MW	104kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using MAP3K14 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using MAP3K14 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.