

MAP4K1 Mouse mAb

Catalog No: #63906

Package Size: #63906-1 50ul #63906-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MAP4K1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Monkey,Rat
Immunogen Description	Purified recombinant fragment of human MAP4K1 (AA:278-477) expressed in mammalian.
Target Name	MAP4K1
Other Names	HPK1
Accession No.	Q92918
Calculated MW	91.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using MAP4K1 mouse mAb against Jurkat (1), A549 (2), mouse kidney (3),COS7 (4),C6 (5) and NIH/3T3 (6) cell lysate.



Flow cytometric analysis of Jurkat cells using MAP4K1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using MAP4K1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using MAP4K1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.