

MAPK10 Mouse mAb

Catalog No: #63908

Package Size: #63908-1 50ul #63908-2 100ul


Description

Product Name	MAPK10 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human MAPK10 (AA: 180-329) expressed in E. Coli.
Target Name	MAPK10
Other Names	JNK3; JNK3A; PRKM10; SAPK1b; p493F12; p54bSAPK
Accession No.	P53779
Calculated MW	52.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using MAPK10 mouse mAb against HEK293 (1), Hela (2), SK-N-SH (3), MCF-7 (4), Jurkat (5), and C6 (6) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.