

MAPK7 Mouse mAb

Catalog No: #63909

Package Size: #63909-1 50ul #63909-2 100ul


Description

Product Name	MAPK7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MAPK7 (AA: 411-556) expressed in E. Coli.
Target Name	MAPK7
Other Names	BMK1; ERK4; ERK5; PRKM7
Accession No.	Q13164
Calculated MW	88.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using MAPK7 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using MAPK7 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.