

MAPK9 Mouse mAb

Catalog No: #63911

Package Size: #63911-1 50ul #63911-2 100ul


Description

Product Name	MAPK9 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MAPK9 (AA: 227-382) expressed in E. Coli.
Target Name	MAPK9
Other Names	JNK2; SAPK; p54a; JNK2A; JNK2B; PRKM9; JNK-55; SAPK1a; JNK2BETA; p54aSAPK; JNK2ALPHA
Accession No.	P45984
Calculated MW	48.1kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using MAPK9 mouse mAb against HEK293 (1) and U251 (2) cell lysate.



Flow cytometric analysis of HeLa cells using MAPK9 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.