

MATN1 Mouse mAb

Catalog No: #63914

Package Size: #63914-1 50ul #63914-2 100ul


Description

Product Name	MATN1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MATN1 expressed in E. Coli.
Target Name	MATN1
Other Names	MAT; SAMS; MATA1; SAMS1; MAT1A
Accession No.	P21941
Calculated MW	54kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Immunofluorescence analysis of HeLa cells using MATN1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.