

MCL-1 Mouse mAb

Catalog No: #63918

Package Size: #63918-1 50ul #63918-2 100ul


Description

Product Name	MCL-1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MCL-1 expressed in E. Coli.
Target Name	MCL-1
Other Names	EAT, MCL1L, MCL1S
Accession No.	Q07820
Calculated MW	37kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using MCL1 mouse mAb against Hela (1), BCBL-1 (2), Jurkat (3) and HL60 (4) cell lysate.



Confocal Immunofluorescence analysis of HepG2 cells using MCL1 mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human lymphnode tissues using MCL1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.