

MCM3 Mouse mAb

Catalog No: #63921

Package Size: #63921-1 50ul #63921-2 100ul


Description

Product Name	MCM3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human MCM3 expressed in E. Coli.
Target Name	MCM3
Other Names	HCC5; P1.h; RLFB; P1-MCM3
Accession No.	P25205
Calculated MW	91kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/50 - 1/200FC:1/200 - 1/400


Images




Western blot analysis using MCM3 mouse mAb against PC-12 (1), HepG2 (2),HeLa (3) and COS7 (4) cell lysate.




Immunofluorescence analysis of HeLa cells using MCM3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using MCM3 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using MCM3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MCM3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.