

MDM4 Mouse mAb

Catalog No: #63922

Package Size: #63922-1 50ul #63922-2 100ul


Description

Product Name	MDM4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MDM4 expressed in E. Coli.
Target Name	MDM4
Other Names	HDMX; MDMX; MRP1; MDM4
Accession No.	O15151
Calculated MW	55kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 ELISA:1/10000


Images




Western blot analysis using MDM4 mouse mAb against HeLa (1), A549 (2) and A431 (3) cell lysate.



Confocal Immunofluorescence analysis of HeLa (left) and L-02 (right) cells using MDM4 mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin.



Immunohistochemical analysis of paraffin-embedded human cerebra (left) and lung carcinoma (right) tissues, showing nuclear localization with DAB staining using MDM4 mouse mAb.



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using MDM4 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.