

MEN1 Mouse mAb

Catalog No: #63927

Package Size: #63927-1 50ul #63927-2 100ul


Description

Product Name	MEN1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MEN1 (AA: 392-554) expressed in E. Coli.
Target Name	MEN1
Other Names	MEAI; SCG2
Accession No.	O00255
Calculated MW	68kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using MEN1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MEN1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.