

MLL Mouse mAb

Catalog No: #63935

Package Size: #63935-1 50ul #63935-2 100ul


Description

Product Name	MLL Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of MLL (aa3751-3968) expressed in E. Coli.
Target Name	MLL
Other Names	MLL
Accession No.	Q03164
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung cancer (left) and esophagus cancer (right), showing nuclear weak staining with DAB staining using MLL mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.