

MLXIPL Mouse mAb

Catalog No: #63936

Package Size: #63936-1 50ul #63936-2 100ul


Description

Product Name	MLXIPL Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MLXIPL (AA: 18-143) expressed in E. Coli.
Target Name	MLXIPL
Other Names	MIO; CHREBP; MONDOB; WBSR14; WS-bHLH; bHLHd14
Accession No.	Q9NP71
Calculated MW	93.1kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using MLXIPL mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.