

## MMP3 Mouse mAb

Catalog No: #63938

Package Size: #63938-1 50ul #63938-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | MMP3 Mouse mAb                                                                           |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | FC                                                                                       |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human MMP3 expressed in E. Coli.                        |
| Target Name           | MMP3                                                                                     |
| Other Names           | SL-1; STMY; STR1; CHDS6; MMP-3; STMY1                                                    |
| Accession No.         | P08254                                                                                   |
| Calculated MW         | 54kDa                                                                                    |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide                                         |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

FC:1/200 - 1/400ELISA:1/10000

## Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of NIH/3T3 cells using MMP3 mouse mAb (green) and negative control (red).

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.