

MOB1A Mouse mAb

Catalog No: #63939

Package Size: #63939-1 50ul #63939-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MOB1A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MOB1A (AA: 1-216) expressed in E. Coli.
Target Name	MOB1A
Other Names	MOB1; MATS1; Mob4B; C2orf6; MOBK1B; MOBKL1B
Accession No.	Q9H8S9
Calculated MW	25kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using MOB1A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.