

mouse Lplunc1 Mouse mAb

Catalog No: #63941

Package Size: #63941-1 50ul #63941-2 100ul


Description

Product Name	mouse Lplunc1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Purified recombinant fragment of mouse Lplunc1 (AA: 248-475) expressed in E. Coli.
Target Name	mouse Lplunc1
Other Names	Bpifb1; RP23-154J12.1
Accession No.	Q61114
Calculated MW	52.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using Lplunc1 mouse mAb against NIH3T3 (1) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.