

Mouse TUG Mouse mAb

Catalog No: #63943

Package Size: #63943-1 50ul #63943-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Mouse TUG Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Purified recombinant fragment of TUG expressed in E. Coli.
Target Name	Mouse TUG
Other Names	ASPCR1, ASPL, ASPS, RCC17, TUG, UBXD9, UBXN9
Accession No.	Q9BZE9
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using TUG mouse mAb against NIH/3T3 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.