

MPS1 Mouse mAb

Catalog No: #63944

Package Size: #63944-1 50ul #63944-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MPS1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of MPS1 expressed in E. Coli.
Target Name	MPS1
Other Names	RPS27
Accession No.	P42677
Calculated MW	95kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Immunohistochemical analysis of paraffin-embedded human lung cancer (left) and human brain (right) tissues using MPS1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.