

MSI2 Mouse mAb

Catalog No: #63953

Package Size: #63953-1 50ul #63953-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MSI2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MSI2 expressed in E. Coli.
Target Name	MSI2
Other Names	MSI2H; MGC3245; FLJ36569; MSI2
Accession No.	Q96DH6
Calculated MW	35kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using MSI2 mouse mAb against NTERA-2 (1), SW620 (2) and T47D (3) cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer (left) and cerebellum tissues (right) using MS12 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.