

MSLN Mouse mAb

Catalog No: #63954

Package Size: #63954-1 50ul #63954-2 100ul


Description

Product Name	MSLN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human MSLN (AA: 37-286) expressed in E. Coli.
Target Name	MSLN
Other Names	MPF; SMRP
Accession No.	Q13421
Calculated MW	69kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using MSLN mouse mAb against CHO3D10 (1), SK-MES-1 (2), and HL-60 (3) cell lysate.



Flow cytometric analysis of HeLa cells using MSLN mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.