

MSN Mouse mAb

Catalog No: #63956

Package Size: #63956-1 50ul #63956-2 100ul


Description

Product Name	MSN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human MSN expressed in E. Coli.
Target Name	MSN
Other Names	moesin
Accession No.	P26038
Calculated MW	67.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using MSN mouse mAb against HeLa (1), A431 (2),Jurkat(3), HEK293(4), and COS7 (5) cell lysate.



Flow cytometric analysis of Jurkat cells using MSN mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon tissues using MSN mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.