

MTHFR Mouse mAb

Catalog No: #63961

Package Size: #63961-1 50ul #63961-2 100ul


Description

Product Name	MTHFR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human MTHFR expressed in E. Coli.
Target Name	MTHFR
Other Names	MTHFR
Accession No.	P42898
Calculated MW	75kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using MTHFR mouse mAb against Rat Heart cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung cancer using MTHFR mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human intima cancer tissues (left) and prostate tissues (right) using MTHFR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.