

MUC12 Mouse mAb

Catalog No: #63964

Package Size: #63964-1 50ul #63964-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MUC12 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MUC12 (AA: extra 371-592) expressed in HEK293-6e cells supernatant.
Target Name	MUC12
Other Names	MUC11; MUC-11;
Accession No.	Q9UKN1
Calculated MW	55.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000FC: 1/200 - 1/400ELISA: 1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using MUC12 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using MUC12 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MUC12 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.