

MUC16 Mouse mAb

Catalog No: #63967

Package Size: #63967-1 50ul #63967-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MUC16 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MUC16 (AA: extra 1-200) expressed in E. Coli.
Target Name	MUC16
Other Names	CA125
Accession No.	Q8WXI7
Calculated MW	1519kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of SK-OV-3 cells using MUC16 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using MUC16 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MUC16 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.