

MUM1 Mouse mAb

Catalog No: #63971

Package Size: #63971-1 50ul #63971-2 100ul


Description

Product Name	MUM1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human MUM1 expressed in E. Coli.
Target Name	MUM1
Other Names	MUM-1; EXPAND1; HSPC211; FLJ14868; FLJ22283; MGC131891; MGC163315; MUM1
Accession No.	Q2TAK8
Calculated MW	78.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Immunohistochemical analysis of paraffin-embedded human melanoma (A), non-Hodgkin's lymphoma (B), using MUM1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human Uterus tissues using MUM1 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.