

MuRF1 Mouse mAb

Catalog No: #63972

Package Size: #63972-1 50ul #63972-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MuRF1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Synthesized peptide of human MuRF1 (AA: 293-304).
Target Name	MuRF1
Other Names	TRIM63; IRF; SMRZ; MURF2; RNF28
Accession No.	Q969Q1
Calculated MW	40.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 FC:1/200 - 1/400 ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HeLa cells using MuRF1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Hela cells using MuRF1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using MuRF1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human muscle tissues using MuRF1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.