

MYL2 Mouse mAb

Catalog No: #63975

Package Size: #63975-1 50ul #63975-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MYL2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of MYL2 expressed in E. Coli.
Target Name	MYL2
Other Names	MLC2; CMH10; DKFZp779C0562
Accession No.	P10916
Calculated MW	19kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using MYL3 (1) and MYL2 (2) mouse mAb against rat fetal heart tissue lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.