

MYL3 Mouse mAb

Catalog No: #63976

Package Size: #63976-1 50ul #63976-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYL3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of MYL3 expressed in E. Coli.
Target Name	MYL3
Other Names	CMH8; VLC1; MLC1V; MLC1SB
Accession No.	P08590
Calculated MW	22kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using MYL3 (1) and MYL2 (2) mouse mAb against rat fetal heart tissues lysate.



Immunohistochemical analysis of paraffin-embedded human heart tissues using MYL3 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.