

NANOG Mouse mAb

Catalog No: #63982

Package Size: #63982-1 50ul #63982-2 100ul


Description

Product Name	NANOG Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of NANOG (aa20-166) expressed in E. Coli.
Target Name	NANOG
Other Names	NANOG
Accession No.	Q9H9S0
Calculated MW	35kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using NANOG mouse mAb against NTERA-2 cell lysate (2).



Confocal Immunofluorescence analysis of NTERA-2 cells (left) and HeLa cells (right) using Nanog mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.