

NAPSA Mouse mAb

Catalog No: #63983

Package Size: #63983-1 50ul #63983-2 100ul


Description

Product Name	NAPSA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human NAPSA (AA: 20-158) expressed in E. Coli.
Target Name	NAPSA
Other Names	KAP; Kdap; NAP1; NAPA; SNAPA
Accession No.	O96009
Calculated MW	45.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using NAPSA mouse mAb against rat liver tissue lysate.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using NAPSA mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using NAPSA mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.