

NCAM1 Mouse mAb

Catalog No: #63985

Package Size: #63985-1 50ul #63985-2 100ul


Description

Product Name	NCAM1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human *** (AA: 568-708) expressed in E. Coli.
Target Name	NCAM1
Other Names	CD56; NCAM; MSK39
Accession No.	P13591
Calculated MW	94.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using NCAM1 mouse mAb against NIH/3T3 (1) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.