

NCOR1 Mouse mAb

Catalog No: #63988

Package Size: #63988-1 50ul #63988-2 100ul


Description

Product Name	NCOR1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of NCOR1 (aa1-192) expressed in E. Coli.
Target Name	NCOR1
Other Names	N-CoR; TRAC1; hN-CoR
Accession No.	O75376
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human cerebra (left) and breast carcinoma tissue (right), showing nuclear location with DAB staining using NCOR1 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.