

NEDD8 Mouse mAb

Catalog No: #63990

Package Size: #63990-1 50ul #63990-2 100ul


Description

Product Name	NEDD8 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NEDD8 expressed in E. Coli.
Target Name	NEDD8
Other Names	NEDD-8
Accession No.	Q15843
Calculated MW	9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using NEDD8 mouse mAb against C6 (1) and HeLa (2) cell lysate.



Immunofluorescence analysis of HeLa cells using NEDD8 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of HeLa cells using NEDD8 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using NEDD8 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.