

NEDD9 Mouse mAb

Catalog No: #63991

Package Size: #63991-1 50ul #63991-2 100ul


Description

Product Name	NEDD9 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human NEDD9 (AA: 82-398) expressed in E. Coli.
Target Name	NEDD9
Other Names	CAS2; CASL; HEF1; CAS-L; CASS2
Accession No.	Q14511
Calculated MW	92.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using NEDD9 mouse mAb against MCF-7 (1), HeLa (2), C2C12 (3),and Hek293 (4) cell lysate.




Immunofluorescence analysis of HeLa cells using NEDD9 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HepG2 cells using NEDD9 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using NEDD9 mouse mAb (green) and negative control (red).



Flow cytometric analysis of K562 cells using NEDD9 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.