

NEFM Mouse mAb

Catalog No: #63994

Package Size: #63994-1 50ul #63994-2 100ul


Description

Product Name	NEFM Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NEFM (AA: 381-443) expressed in E. Coli.
Target Name	NEFM
Other Names	NFM; NEF3; NF-M
Accession No.	P07197
Calculated MW	102.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using NEFM mouse mAb against NTERA-2 (1), SK-N-SH (2), and PC-12 (3) cell lysate.



Flow cytometric analysis of SK-N-SH cells using NEFM mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.