

NEFM Mouse mAb

Catalog No: #63995

Package Size: #63995-1 50ul #63995-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NEFM Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NEFM (AA: 779-916) expressed in E. Coli.
Target Name	NEFM
Other Names	NFM; NEF3; NF-M
Accession No.	P07197
Calculated MW	102.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/50 - 1/200FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using NEFM mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Hela cells using NEFM mouse mAb (green) and negative control (red).



Flow cytometric analysis of Raji cells using NEFM mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.