

Neurod Mouse mAb

Catalog No: #63996

Package Size: #63996-1 50ul #63996-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Neurod Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human Neurod (AA: 26-91) expressed in E. Coli.
Target Name	Neurod
Other Names	NEUROD1; BETA2; BHF-1; MODY6; bHLHa3
Accession No.	Q13562
Calculated MW	40kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Western blot analysis using Neurod mouse mAb against WERI-Rb-1 cell lysate.



Flow cytometric analysis of Hela cells using Neurod mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human brain tissues using Neurod mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using Neurod mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.