

NFKB2 Mouse mAb

Catalog No: #63998

Package Size: #63998-1 50ul #63998-2 100ul


Description

Product Name	NFKB2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human NFKB2 (AA: 712-900) expressed in E. Coli.
Target Name	NFKB2
Other Names	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-kB2; p49/p100
Accession No.	Q00653
Calculated MW	96.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using NFKB2 mouse mAb against NIH/3T3 (1), MCF-7 (2), and THP-1 (3) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.