

NGFR Mouse mAb

Catalog No: #64002

Package Size: #64002-1 50ul #64002-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NGFR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NGFR (AA: 1-200) expressed in E. Coli.
Target Name	NGFR
Other Names	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR
Accession No.	P08138
Calculated MW	45.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using NGFR mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebrum tissues using NGFR mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using NGFR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.