

NME1 Mouse mAb

Catalog No: #64004

Package Size: #64004-1 50ul #64004-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NME1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NME1 expressed in E. Coli.
Target Name	NME1
Other Names	NB; AWD; NBS; GAAD; NM23; NDPKA; NDPK-A; NM23-H1; NME1
Accession No.	P15531
Calculated MW	17kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Immunofluorescence analysis of HeLa cells using NME1 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of Jurkat cells using NME1 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Placenta tissues using anti-NME1 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.