

NPM2 Mouse mAb

Catalog No: #64013

Package Size: #64013-1 50ul #64013-2 100ul


Description

Product Name	NPM2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human NPM2 (AA: 1-214) expressed in E. Coli.
Target Name	NPM2
Accession No.	Q86SE8
Calculated MW	24.1kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/100 - 1/400 ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using NPM2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissues using NPM2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human brain tissues using NPM2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.