

NT5E Mouse mAb

Catalog No: #64026

Package Size: #64026-1 50ul #64026-2 100ul


Description

Product Name	NT5E Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Rat,Rabbit
Immunogen Description	Purified recombinant fragment of human NT5E (AA: 30-250) expressed in E. Coli.
Target Name	NT5E
Other Names	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA
Accession No.	P21589
Calculated MW	63.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using NT5E mouse mAb against A431 cell lysate.




Flow cytometric analysis of A431 cells using NT5E mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using NT5E mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using NT5E mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded Mouse testis using NT5E mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat kidney using NT5E mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit testis using NT5E mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit spleen using NT5E mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.