

NTRK3 Mouse mAb

Catalog No: #64028

Package Size: #64028-1 50ul #64028-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NTRK3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant extracellular fragment of human NTRK3 (aa32-429) fused with hlgGfc tag expressed in HEK293 cells.
Target Name	NTRK3
Other Names	TRKC; gp145(trkC)
Accession No.	Q16288
Calculated MW	94kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000ELISA: 1/10000

Images



Western blot analysis using NTRK3 mouse mAb against extracellular domain of human NTRK3 (aa32-429).



Immunohistochemical analysis of paraffin-embedded human lymph node using NTRK3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.