

OLIG2 Mouse mAb

Catalog No: #64029

Package Size: #64029-1 50ul #64029-2 100ul


Description

Product Name	OLIG2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human OLIG2 expressed in E. Coli.
Target Name	OLIG2
Other Names	BHLHB1; OLIGO2; RACK17; PRKCBP2; bHLHe19; OLIG2
Accession No.	Q13516
Calculated MW	32kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000 ICC:1/200 - 1/1000 ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunofluorescence analysis of U251 cells using OLIG2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human muscle tissues using OLIG2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.