

ONECUT3 Mouse mAb

Catalog No: #64030

Package Size: #64030-1 50ul #64030-2 100ul


Description

Product Name	ONECUT3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ONECUT3 (AA: 350-455) expressed in E. Coli.
Target Name	ONECUT3
Other Names	OC3
Accession No.	O60422
Calculated MW	50kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of A431 cells using ONECUT3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using ONECUT3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ONECUT3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.