

P16 (Mouse and Human) Mouse mAb

Catalog No: #64037

Package Size: #64037-1 50ul #64037-2 100ul


Description

Product Name	P16 (Mouse and Human) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC
Species Reactivity	Human,Mouse,Rat,Rabbit
Immunogen Description	Purified recombinant fragment of P16 expressed in E. Coli.
Target Name	P16 (Mouse and Human)
Other Names	P16
Accession No.	P42771
Calculated MW	16.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images




Immunohistochemical analysis of paraffin-embedded rat liver tissue (A), human brain tissue (B) and brain tumor (C), showing nuclear localization using P16 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Mouse liver(A)Mouse testis(B)Mouse kidney(C) using P16 (Mouse and Human) mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat liver(A)Rat testis(B)Rat kidney(C) using P16 (Mouse and Human) mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit liver(A)Rabbit testis(B) using P16 (Mouse and Human) mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.