

P2RY1 Mouse mAb

Catalog No: #64039

Package Size: #64039-1 50ul #64039-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	P2RY1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IHC;FC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human P2RY1 (AA: extra mix) expressed in E. Coli.
Target Name	P2RY1
Other Names	P2Y1
Accession No.	P47900
Calculated MW	42kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using P2RY1 mouse mAb against SPC-A-1 (1) and C6 (2) cell lysate.




Flow cytometric analysis of HeLa cells using P2RY1 mouse mAb (green) and negative control (red).



Flow cytometric analysis of K562 cells using P2RY1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using P2RY1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using P2RY1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.