

P2RY12 Mouse mAb

Catalog No: #64040

Package Size: #64040-1 50ul #64040-2 100ul


Description

Product Name	P2RY12 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human P2RY12 (AA: extra mix) expressed in E. Coli.
Target Name	P2RY12
Other Names	HORK3; P2Y12; ADPG-R; BDPLT8; SP1999; P2T(AC); P2Y(AC); P2Y(12)R; P2Y(ADP); P2Y(cyc)
Accession No.	Q9H244
Calculated MW	39.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using P2RY12 mouse mAb against PC-3 (1) and C6 (2) cell lysate.



Flow cytometric analysis of Hela cells using P2RY12 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.