

P2RY2 Mouse mAb

Catalog No: #64042

Package Size: #64042-1 50ul #64042-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	P2RY2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human P2RY2 (AA: extra mix) expressed in E. Coli.
Target Name	P2RY2
Other Names	P2U; HP2U; P2U1; P2UR; P2Y2; P2RU1; P2Y2R
Accession No.	P41231
Calculated MW	42.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using P2RY2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.